

Special examinations of milk- and other powders

Solumixer

To determine the solubility of various powders.
Complies with the test specifications of ADPI and DLG.
Stirring blades made of stainless steel.
Complete with a blender glass.



09.0040

Order No.	Type	Qty per unit
09.0040-230CH	Solumixer Universal 230 VAC/50 Hz with Swiss power plug digital display, one fixed program (3600 rpm, 90 sec) for solubility and variable programming for general use (1800-3600 rpm, timer 1 - 9 min)	1
09.0040-230CEE	Solumixer Universal 230 VAC/50 Hz with Schuko power plug digital display, one fixed program (3600 rpm, 90 sec) for solubility and variable programming for general use (1800-3600 rpm, timer 1 - 9 min)	1
09.0040-115	Solumixer Universal 115 VAC/60 Hz with US power plug digital display, one fixed program (3600 rpm, 90 sec) for solubility and variable programming for general use (1800-3600 rpm, timer 1 - 9 min)	1
09.0041-230	Solumixer Standard 230 VAC/50 Hz, without digital display, with a fixed program (3600 rpm, 90 sec) for solubility	1
09.0041-115	Solumixer Standard 115 VAC/60 Hz without digital display, with a fixed program (3600 rpm, 90 sec) for solubility	1

Accessories for Solumixer

Order No.	Type	Qty per unit
09.0051	Spare blender glass	1
09.0052	Spare stirrer blade	1
09.0053	Spare driving belt	1
09.0056	Spare rocker arm	1

Waring Mixer

Robust and efficient laboratory mixer for food samples.
Blender glass: 1.0 liters
Stand with knife block

Specifications:

Speed : 500 - 2000 rpm
Power: 840 W
Voltage : 230 VAC

Order No.	Type	Qty per unit
35.0100	Waring mixer LB 20 EG complete with blender glass	1
35.0122	Spare blender glass	1