

# TERRARACK TIPS: THE MOST ECO-FRIENDLY PIPETTE TIP RACK AROUND



ONE TREE PLANTED FOR EVERY 20 BOXES YOU PURCHASE.

### LESS WEIGHT, LESS WASTE, 100% RECYCLABLE

TerraRack is a revolutionary new concept in tip racks. Sturdy as conventional racks, yet made with less than half the plastic and completely recyclable. In fact, the hinged TerraRack shell is made from PETE, which is easily recycled and in high demand.

Single-use and non-refillable, TerraRacks can be used independently or with a reusable TerraBase for added support and stability.

For complete convenience and cleanliness, all TerraRacks are pre-sterilized, eliminating the need to autoclave (TerraRacks are not autoclavable).

TerraRack is available in LTS and universal-fit tips.

### EVERY RACK IS FRESH, STERILE

TerraRacks are pre-sterilized. Unlike refill systems, where there is a potential for placing clean tips into a contaminated rack, with TerraRack you can be certain that you're working with a fresh, clean rack every time.

### 100% RECYCLABLE

In addition to being completely recyclable, TerraRack is largely composed of PETE, which is easily recycled and aggressively recovered from the plastic waste stream. Recycled PETE is in high demand, particularly by the beverage industry.

### LESS WASTE

In addition to containing much less plastic than conventional racks, TerraRacks are extremely compressible, which means that used racks occupy significantly less space in a waste container or storage vessel.

### LESS WEIGHT

TerraRacks weigh less than half that of conventional injection-molded tip racks. As such, TerraRacks cost significantly less to transport and, where disposal rates are based on weight, cost less to throw away.

		Nominal Pipette Volumes								
LTS Models		2 µL	10 µL	20 µL	50 µL	100 µL	200 µL	300 µL	1000 µL	1200 µL
	LTS									
	LTS Pre-sterilized	●	●	●	●	●	●	●	●	●
	LTS Pre-sterilized, filter	●	●	●	●	●	●	●	●	●
Universal-fit Models		2 µL	10 µL	20 µL	25 µL	100 µL	200 µL	300 µL	1000 µL	1200 µL
	UNV									
	UNV Pre-sterilized	● <sup>S</sup>	● <sup>S</sup>	●	●	●	●	●	●	●
	UNV Pre-sterilized, filter	● <sup>S</sup>	● <sup>S</sup>	●	●	●	●	●	●	●

<sup>S</sup> Available in ShiftGard

CODE	DESCRIPTION	PRICE
<b>TERRARACK WITH LTS TIPS – FOR RAININ LTS PIPETTES</b>		
17014960	TR-L10S (960 tips in 10 racks), Presterile 20 µL	\$ 98.00
17014961	TR-L10F (960 tips in 10 racks), Filter 20 µL	\$ 110.00
17014962	TR-L250S (960 tips in 10 racks), Presterile 200 µL	\$ 98.00
17014963	TR-L200F (960 tips in 10 racks), Filter 200 µL	\$ 110.00
17014964	TR-L300S (768 tips in 8 racks), Filter 300 µL	\$ 98.00
17014965	TR-L300F (768 tips in 8 racks), Filter 300 µL	\$ 110.00
17014966	TR-L1000S (768 tips in 8 racks), Presterile 1000 µL	\$ 98.00
17014967	TR-L1000F (768 tips in 8 racks), Filter 1000 µL	\$ 110.00
17014968	TR-L1200S (768 tips in 8 racks), Presterile 1200 µL	\$ 98.00
17014969	TR-L1200F (768 tips in 8 racks), Filter 1200 µL	\$ 110.00
<b>TERRARACK WITH TRADITIONAL TIPS – FOR UNIVERSAL PIPETTES</b>		
17014970	TR-10S (960 tips in 10 racks), Presterile 10 µL	\$ 98.00
17014971	TR-10F (960 tips in 10 racks), Filter 10 µL	\$ 110.00
17014972	TR-10GS (960 tips in 10 racks), ShaftGard Presterile 10 µL	\$ 98.00
17014973	TR-10GF (960 tips in 10 racks), ShaftGard Filter 10 µL	\$ 110.00
17014974	TR-20F (960 tips in 10 racks), Filter 30 µL	\$ 110.00
17014975	TR-100F (960 tips in 10 racks), Filter 100 µL	\$ 110.00
17014976	TR-250S (960 tips in 10 racks), Presterile 250 µL	\$ 98.00
17014977	TR-200F (960 tips in 10 racks), Filter 200 µL	\$ 110.00
17014978	TR-300S (768 tips in 8 racks), Presterile 300 µL	\$ 98.00
17014979	TR-300F (768 tips in 8 racks), Filter 300 µL	\$ 110.00
17014980	TR-1000S (768 tips in 8 racks), Presterile 1000 µL	\$ 98.00
17014981	TR-1000F (768 tips in 8 racks), Filter 1000 µL	\$ 110.00
<b>OPTIONAL TERRA BASE – REUSABLE BASE FOR TERRARACK. ONLY REQUIRED IF ADDITIONAL STABILITY IS DESIRED.</b>		
<b>Non-slip accessory.</b>		
• Each TerraBase provides secure footing on slippery bench surfaces for an enhanced pipetting experience.		
<b>Durable construction.</b>		
• Polycarbonate composition provides high impact resistance. TerraBases are autoclavable and chemical resistant.		
17014982	TR-SB TerraRack Base - short, 1 base to fit 10-300µL TerraRack racks	\$ 37.70
17014983	TR-TB TerraRack Base - tall, 1 base to fit 1000-1200µL TerraRack racks	\$ 41.30