

Only \$1,742
Save \$436!

Rainin NanoRep

World's Versatile Electronic Repeater Pipette

With the new Rainin NanoRep, life science researchers can accurately dispense solutions from 50 mL to as low as 100 nL, meeting the growing demand for protocol miniaturization and optimization. Researchers can dispense as many as 1,000 aliquots from any NanoRep syringe tip. NanoRep is a practical alternative to expensive benchtop systems, combining the versatility and value of a repeater pipette with a unique non-contact dispense option that ensures repeatable accuracy, improves efficiency and reduces the risk of cross-contamination.

Viscous, volatile, and other solutions can be pipetted accurately and safely due to positive displacement.

Rainin NanoRep...The most versatile pipette in your lab!



Request a **FREE QUOTE** or **TRIAL NOW!**
Email onelab@onelab.co.nz or call **0800 ONELAB**
www.mt.com/pipettes

*Valid until 31st of October 2023

RAININ
Pipetting 360°