


# QUALITY MATTERS: FILTER TIPS – BUY 3 PACKS GET ONE FREE\*



### Low-Retention

All Vertex tips are designed to be naturally low-binding, but tips with NoStick® resin are considerably more hydrophobic than standard tips, enabling even viscous fluids to be dispensed more completely.

### Clean Liquid Release

Vertex tips deliver a much cleaner liquid release, dispensing all of the sample and improving pipetting accuracy. Users will maintain complete fluid control, even when working with smaller volumes, while not compromising the accuracy and precision of their pipetting.

**\*DEAL PRICE EQUALS  
\$58.50 PER 10 RACKS!!!**

### Key Features:

Vertex® Pipette Tips were meticulously designed using a specially formulated NoStick® resin to provide researchers with the highest level of liquid handling performance. Vertex Pipette Tips were created with the belief that when it comes to liquid handling, less is more.

Vertex® Pipette Tips are produced under rigorous controls and safeguards- ISO9001 certified Quality Assurance process and state of the art fully automated tip manufacturing, filtering and inspecting is incomparable in our industry.

- Less fluid retention
- Less sample loss
- Less time spent
- Less plastic wasted



**IN STOCK NOW**

DESCRIPTION	CODE	PACK	PRICE – PLACE YOUR ORDER WITH THIS PRICE	BUY 3 PACKS GET ONE FREE MAKES THE PRICE PER PACK
Tip 10 µL Racked, Filter Tips, NoStick Clear, Graduated	SSIB4117NSF	10 Racks	\$ 78.00	<b>\$ 58.50</b>
Tip 10 µL XL Racked Filter, NoStick Clear	SSIB4137NSF	10 Racks	\$ 78.00	<b>\$ 58.50</b>
Tip 20 µL Racked Filter NoStick Clear, Graduated	SSIB4237NAF	10 Racks	\$ 78.00	<b>\$ 58.50</b>
Tip 100 µL Racked Filter NoStick Clear, Graduated	SSIB4237NBF	10 Racks	\$ 78.00	<b>\$ 58.50</b>
Tip 200 µL Racked Filter NoStick Clear, Graduated	SSIB4237NSF	10 Racks	\$ 78.00	<b>\$ 58.50</b>
Tip 300 µL Racked Filter, Sterile NoStick Clear, Graduated	SSIB4927NSF	10 Racks	\$ 78.00	<b>\$ 58.50</b>
Tip 1000 µL Racked Filter NoStick Clear, Graduated	SSIB4337NSF	10 Racks	\$ 84.00	<b>\$ 63.00</b>
Tip 1250 µL Racked Filter Sterile NoStick Clear, Graduated	SSIB4347NSF	10 Racks	\$ 88.75	<b>\$ 66.60</b>

**PROMO CODE: "TIP PROMO" must be included on your order.** Mix and match is possible, lowest priced tip pack will be Free in this case. Promo code cannot be applied retrospectively.