


results matter.



Omni Tip Homogenizing Kits

**15% off and Extended 3 year
Warranty on OMNI TH PCR kits.**

Valid until 31st October 2020

OMNI Tip™ Homogenizing Kits come complete with a redesigned TH motor, OMNI Tip™ adapter, 12 hard or soft tissue patented OMNI Tip™ Plastic Generator Probes, and tool kit.

The kits can be used to process samples with a range of <200 µL to 100 mL. The perfect solution for a small laboratory, this kit is cost effective and will handle a wide range of applications for multiple users. Kits are also available with optional stainless steel generator probes.



Features

- Economical: <1% of the cost of stainless steel probes
- Eliminates Cross Contamination
- Ideal for PCR, RNA, and DNA sample preparation
- Maintenance Free

Proprietary OMNI Tip™ probes are available in both soft tissue (deaggregation) and hard tissue versions. Engineered from durable plastics, OMNI Tips™ are reusable (by chemical sterilization or autoclave) and are economical enough to be disposed of after each use to eliminate cross contamination.

The outer tube of OMNI Tips™ is made of polycarbonate. The inner shafts are made of Ultem 1000, a polyetherimide. The OMNI Tips™ patented design eliminates the hassle of disassembly and maintenance of traditional stainless steel generators.

Treatable volume - 250 µL to 30 mL

Page 2 for further details and pricing.

Omni Tip Homogenizing Kits



15% off and Extended 3 year Warranty on OMNI TH PCR kits.

Valid until 31st October 2020

Catalogue Number	Description	List Price	Promo Price
OMNITH220-PCRD	Omni Tip Kit with 220 volt Omni TH & 12 Soft Tissue Omni Tips	\$ 1,935	\$ 1,644
OMNITH220-PCR5D	Omni Tip Kit with 220 volt Omni TH, 12 Soft Tissue Omni Tips & 5x75mm Stainless Steel Generator Probe (G5-75)	\$ 3,274	\$ 2,783
OMNITH220-PCR7D	Omni Tip Kit with 220 volt Omni TH, 12 Soft Tissue Omni Tips & 7x95mm Stainless Steel Generator Probe (G7-115ST)	\$ 3,274	\$ 2,783
OMNITH220-PCR10D	Omni Tip Kit with 220 volt Omni TH, 12 Soft Tissue Omni Tips & 10x95mm Stainless Steel Generator Probe (G10-115)	\$ 3,472	\$ 2,951
OMNITH220-PCRH	Omni Tip Kit with 220 volt Omni TH & 12 Hard Tissue Omni Tips	\$ 1,935	\$ 1,644
OMNITH220-PCR5H	Omni Tip Kit with 220 volt Omni TH, 12 Hard Tissue Omni Tips & 5x75mm Stainless Steel Generator Probe (G5-75)	\$ 3,274	\$ 2,783
OMNITH220-PCR7H	Omni Tip Kit with 220 volt Omni TH, 12 Hard Tissue Omni Tips & 7x95mm Stainless Steel Generator Probe (G7-115ST)	\$ 3,274	\$ 2,783
OMNITH220-PCR10H	Omni Tip Kit with 220 volt Omni TH, 12 Hard Tissue Omni Tips & 10x95mm Stainless Steel Generator Probe (G10-115)	\$ 3,472	\$ 2,951

Prices exclude GST and local freight.

CONTACT US

09 414 5599

onelab@onelab.co.nz