

PIPETTE MANAGEMENT SIMPLIFIED



SmartStand and Hang-Up in One!
SmartStand turns into a "Smart Hang-Up" with an accessory fastening kit.



The Rainin SmartStand is an intelligent, easy-to-use pipette asset management system for busy labs and lab managers who want to:

- Reduce out-of-service/specification-related costs
- Proactively manage compliance
- Save time on record keeping and reporting

Pipette Management Simplified SmartStand Fits Your Workflow By reading the RFID chip that's in every Rainin XLS and XLS+ brand pipette, SmartStand instantly displays the compliance status for up to four pipettes – manual or electronic.

Pick up a pipette and SmartStand's color LCD screen switches to a detailed view of that pipette's service and calibration status.

SmartStand works independently or integrates with EasyDirect™ Pipette Asset Management software to create a laboratory-wide pipette asset management and tracking system.

EasyDirect simplifies lab asset management by enabling users to record pipette usage, flag pipettes for service, transfer data between other lab management systems (LIMS) and by recording unique data (such as owner, quick check methods and inventory number) onto each pipette's RFID chip.

SPECIAL PROMOTION:
Special Introductory Pricing on Rainin Smart Stand Lite Touch System (LTS) Pipette Kits and Easy Direct Software:



SmartStand · Flexible, scalable · For manual and electronic XLS and XLS+ pipettes · User defined service/calibration intervals · QuickCheck reminders for periodic performance verification · Intuitive graphical interface · 110/220V, 50-60 Hz · Connects via Bluetooth* or USB · Stands alone or use optional kits to mount the head to a shelf or wall

EasyDirect Pipette Asset Management · Connects up to multiple SmartStands to a PC · Tracks and manages all brands/models of pipettes · Reads/writes to RFID chips on Rainin XLS and XLS+ pipettes · Export to Excel, csv or xml (append, overwrite or delete) · Automatic export scheduler for autonomous operation · Schedule service requests · System requirements: Windows 7 (SP1) or higher, 5GB of free memory and 4 GB of RAM



CODE	DESCRIPTION	UOM	PRICE
SMARTSTAND PIPETTE KITS – INCLUDES PIPETTES AS INDICATED, WITH SMART STAND!			
30386597	SMARTSTAND PIPETTE KIT L-XLS+ (2. 20. 200.1000 Single Channel LTS Pipettes)	Each	\$ 1,822.22
30386598	SMARTSTAND PIPETTE KIT SL-XLS+ (2. 20. 200.1000 Single Channel LTS Pipettes)	Each	\$ 1,822.22
30386599	SMARTSTAND PIPETTE KIT L-XLS+ (20. 200.1000 Single Channel LTS Pipettes)	Each	\$ 1,600.00
30386730	SMARTSTAND PIPETTE KIT SL-XLS+ (20. 200.1000 Single Channel LTS Pipettes)	Each	\$ 1,600.00
30386733	SMARTSTAND PIPETTE KIT L-XLS+ (20. 200. 1000 Single Channel LTS Pipettes, and an 8-channel L8-200)	Each	\$ 2,693.33
30386734	SMARTSTAND PIPETTE KIT L-XLS+ (20. 200.1000 Single Channel LTS Pipettes, and an 8-channel L8-20)	Each	\$ 2,693.33
30386735	SMARTSTAND PIPETTE KIT L-XLS+ (20. 200.1000 Single Channel LTS Pipettes and a 12-channel L12-200)	Each	\$ 2,835.56
30386736	SMARTSTAND PIPETTE KIT L-XLS+ (20. 200.1000 Single Channel LTS Pipettes and Single Channel LTS Pipettes. And a 12-Channel L12-20)	Each	\$ 2,835.56
30386737	SMARTSTAND PIPETTE KIT L-XLS+ (200.1000. 5000 Single Channel LTS Pipettes)	Each	\$ 1,866.67
30386738	SMARTSTAND PIPETTE KIT L-XLS+ (2.10. 20. 200 Single Channel LTS Pipettes)	Each	\$ 1,822.22
30386739	SMARTSTAND PIPETTE KIT L-XLS+ (Two 8-channel pipettes L8-20 and L8-200)	Each	\$ 2,595.56
30386740	SMARTSTAND PIPETTE KIT L-XLS+ (L12-20 Single Channel LTS Pipettes and a 12-channel L12-200)	Each	\$ 2,924.44
ELECTRONIC PIPETTES			
30386731	SMARTSTAND PIPETTE KIT E4 L-XLS+ (20. 200. 1000 Single Channel Electronic LTS Pipettes)	Each	\$ 3,644.44
30386732	SMARTSTAND PIPETTE KIT SE4 XLS+ (20. 200. 1000 Single Channel Electronic Standard Shaft Pipettes- not LTS, universal fit)	Each	\$ 3,644.44
EASYDIRECT SOFTWARE			
30475915	EasyDirect PAM License - Enterprise	Each	\$5,714.29
30475913	EasyDirect PAM License - Site	Each	\$2,006.08