

# PCR PLATES FOR YOUR INSTRUMENTS



### Key Features:

- Well ridges reduce sample evaporation when sealing with film, foil or associated strip caps
- Ultra-thin, uniform wells ensure optimal heat transfer for high reaction efficiency
- Conveniently packaged in small quantities of 10 plates to minimize risk of contamination
- High contrast black alphanumeric for easy well identification
- Certified free from DNase, RNase, DNA, and PCR inhibitors
- Broad thermal cycler compatibility.
- High-contrast, black alphanumeric indexing for easy well identification.
- Crystal-clear wells to easily verify sample volume.
- Bar-coded semi-skirted and full-skirted plates available.
- Manufactured from high-clarity, ultra-pure polypropylene (meets FDA 21 CFR 177.1520 & USP Class VI).
- For use in medical, pharmaceutical and food industry, and molecular biology.
- Certified free of detectable DNase, RNase, DNA, PCR inhibitors, and tested pyrogen-free.
- Made in the USA

We offer the largest selection of plates that are designed to precisely fit most thermal cyclers on the market. Each is designed to meet and exceed the specifications and performance of plates supplied by the original equipment manufacturer.

**IN STOCK NOW**

DESCRIPTION	CODE	PACK	PRICE
PCR Plate 0.2ml 96-Well semi-skirted (ABI-style) with flat top, Clear	SSIB3420-00	100pcs	\$ 351.40
<b>Compatibility:</b>			
Applied Biosystems® 2700, 3xxx / 3xxx XL / 3xxx DX, 5700, PRISM 7000, 7300, 7500, GeneAmp® 9700			
Bio-Rad® / MJ Research® Gene Cycler			
Eppendorf® Mastercycler ep Gradient			
Biometra® T1 Thermal Cycler, Tgradient, Trobot, Uno, Uno II			
PCR Plate 0.2ml Semi Skirt, Standard 96-well, Clear	SSIB3450-00	100pcs	\$ 351.40
<b>Compatibility:</b>			
Applied Biosystems® 2700, 3xxx / 3xxx XL / 3xxx DX, 5700, PRISM 7000, 7300, 7500, 9600, GeneAmp® 9700			
Bio-Rad® / MJ Research® iCycler, IQ5, MyiQ, PTC-100, PTC-200, PTC-225 Tetra, Dyad			
Eppendorf® Mastercycler, Mastercycler ep Gradient, Mastercycler ep Realplex, Mastercycler Nexus /Gradient, Mastercycler Nexus Eco			
Biometra® T1 Thermal Cycler, Tgradient, Trobot, Uno, Uno II			
PCR Plate 0.2ml 96-well semi-skirted (ABI-style) with top ridge, Clear	SSIB3425-00	100pcs	\$ 351.40
<b>Compatibility:</b>			
Applied Biosystems® 2700, 2720, 6100, 3xxx / 3xxx XL / 3xxx DX, 5700, PRISM 7000, 7300, 7500, 9600, GeneAmp® 9700			
Bio-Rad® / MJ Research® Gene Cycler			
Eppendorf® Mastercycler ep Gradient			

Promotion is valid until 30th August 2020 or while stocks last • Prices exclude GST and local freight