

# STERILISE MEDIA IN MINUTES: MESPRESSO 500 MICRO MEDIAPREPARATOR



## TOTAL TIME OF STERILISATION: 7 MINUTES

- **Mespresso can be installed anywhere.** All you need is an electrical socket.
- **Personal.** Available and ready to use 24/7.
- **Rapid.** Using microwave technology, the melting of the liquid media takes only a short time – as short as it takes to have a cup of coffee – just 7 mins sterilisation time.
- **Touch screen** All functions are operated from the touchscreen – which also displays the progress of the sterilisation process.
- **Uses standard Lab vessels.** Liquid media such as microbiological media, buffer solutions, carbohydrate solutions or water can be sterilised in a standard lab vessels of 100-500ml capacity and temperature resistant up to 135oC (glass or plastic – PMP, PP, PTFE, FEP).
- **Data archiving.** Each sterilisation course is saved on a memory card that can be copied to any computer.



CODE	DESCRIPTION	PRICE	SPECIAL PRICE
SWISPCM500.EU	Mespresso 500, micro Mediaprepator	\$6,995.00	\$6,695.00



## A REVOLUTION IN LIQUID MEDIA STERILISING

All Equipment undergoes electrical safety & function checks at Onelab technical service department prior to dispatch, to ensure maximum safety and satisfaction