

# IKA RCT BASIC – THE BEST SELLING HOTPLATE STIRRER FOR THE LABORATORY



## IKA

designed for scientists

- Strong motor for a higher range of speeds
- Additional temperature control mode for faster heating of medium
- Includes PT 1000 temperature sensor (PT 1000.60)

DESCRIPTION	CODE	PRICE	PROMO PRICE
IKA RCT Basic, Magnetic Stirring with Heating incl. PT1000	IKAA20003581	\$1,699	<b>\$899</b>

Promotion is valid until 30th August 2020 or while stocks last • Prices exclude GST and local freight

### Key Features:

- Integrated temperature control
- Incl. PT 1000 temperature sensor (PT 1000.60)
- Exact temperature and speed setting via digital display, even when switched off
- Set safety temperature limit displayed digitally
- Hot Top indicator > hot surface warning to prevent burns!
- Digital error code display
- With adjustable safety circuit of heating plate temperature (50 - 360 °C)
- Safety magnetic stirrer with heating, suitable for unsupervised operation
- Bushing according to DIN 12878 for connecting a contact thermometer, e.g. ETS-D5, enables precise temperature control
- High level of safety thanks to improved heat control technology
- Enclosed assembly (IP 42) guarantees long service life
- Aluminium heating plate for optimum heat transfer
- Improved magnetic adhesion
- Incl. protection cover H 100
- Max. stirring quantity (H<sub>2</sub>O): 20 l
- Speed range: 50 - 1500 rpm
- Temperature range: RT - 310 °C

All Equipment undergoes electrical safety & function checks at Onelab technical service department prior to dispatch, to ensure maximum safety and satisfaction